
# LAP-BAND<sup>®</sup> System Tubing Kit

DIRECTIONS FOR USE (DFU)

Rx Only

LAP-BAND® System Tubing Kit (English) .....	1
Комплект тръбички за система LAP-BAND® (Български) .....	3
Sada hadiček systému LAP-BAND® (Česky) .....	5
LAP-BAND® -systemets slangesæt (Dansk) .....	7
Slangenset van het LAP-BAND®-systeem (Nederlands) .....	9
LAP-BAND®-i süsteemi torukomplekt (Eesti keel) .....	11
Kit de tubulure du système LAP-BAND® (Français).....	13
LAP-BAND®-System-Füllschlauch-Kit (Deutsch) .....	15
Kit σωληνώσεων για το σύστημα LAP-BAND® (Ελληνικά) .....	17
LAP-BAND® rendszer tömlőkészle (Magyar).....	19
Kit di tubi del sistema LAP-BAND® (Italiano) .....	21
<b>LAP-BAND® 시스템 튜브 키트 (한국어) .....</b>	<b>23</b>
LAP-BAND® Cauruļu komplekts (Latviski).....	25
LAP-BAND® sistemas vamzdelių rinkinys (Lietuviškai).....	27
Slangesett for LAP-BAND®-system (Norsk).....	29
Zestaw rurek systemu LAP-BAND® (Polski).....	31
Kit de Tubo do Sistema LAP-BAND® (Português).....	33
Kitul cu tubulatură pentru sistemul LAP-BAND® (Română).....	35
Комплект трубок системы LAP-BAND® (русский) .....	37
Súprava hadičky pre systém LAP-BAND® (slovenčina) .....	39
Komplet cevki za sistem LAP-BAND® (Slovenščina) .....	41
Kit de tubo del sistema LAP-BAND® (Español) .....	43
LAP-BAND®-systemets slangats (Svenska) .....	45
LAP-BAND® Sistemi Hortum Kiti (Türkçe) .....	47
Трубовий комплект системи LAP-BAND® (українська) .....	49

# LAP-BAND® System Tubing Kit

## INTRODUCTION

The LAP-BAND® Adjustable Gastric Banding System is designed to induce weight loss in severely obese patients by limiting food consumption. With this system, no cutting or stapling of the stomach is required and there is no bypassing of portions of the stomach or intestines.

The LAP-BAND® System Tubing Kit is a replacement part for the LAP-BAND® Adjustable Gastric Banding System. Deflation of the band may occur due to leakage at the Stainless Steel Connector transition or the connecting Silicone Tubing. If such complications should arise, the individual component can be replaced with the LAP-BAND® System Tubing Kit.

LAP-BAND® System Tubing Kit  
(Ref. No. B-20401)

## DESCRIPTION

The LAP-BAND® System Tubing Kit contains the following components:

- 1 each Silicone Tubing, 50.8 cm (sterile)
- 2 each Stainless Steel Connectors (sterile)
- 1 each End Plug (sterile)
- 1 each Band Priming Needle, 40.5 mm (sterile)

## INDICATIONS

Indications for Silicone Tubing or Stainless Steel Connector replacement are:

- Repositioning or relocating of the Access Port
- A tubing leak at or near the Stainless Steel Connector

## CONTRAINDICATIONS

The LAP-BAND® System Tubing Kit is contraindicated in patients where the LAP-BAND® System is contraindicated, and:

1. Patients who have an infection anywhere in the body or where the possibility of contamination prior to or during the surgery exists.
2. Patients who are known to have, or are suspected to have, an allergic reaction to materials contained in the system or who have exhibited a pain intolerance to implanted devices.

## COMPLICATIONS

Complications which may result from the use of this product include the risks associated with the medications and methods utilized in the surgical procedure, the risks associated with any surgical procedure and the patient's degree of intolerance to any foreign object implanted in the body.

Infection can occur in the immediate postoperative period or years after insertion of the device. In the presence of infection or contamination, removal of the device may be indicated. Deflation of the band may occur due to leakage from the band, the port, or the connecting tubing.

The material in this device has been shown in biocompatibility studies to cause slight irritation in intermuscular implantation in animal models.

## RISK ASSOCIATED WITH RE-USE

The LAP-BAND® Tubing Kit is not intended to be re-sterilized or re-used. Cleaning and autoclaving processes may cause damage to tubing and/or components resulting in improper function and possible failure. Reuse of the tubing and/or components can cause infection to the patient.

## HOW SUPPLIED

The LAP-BAND® System Tubing Kit components are sterile and for single use only.

The LAP-BAND® System Tubing Kit components are provided sterile in double packaging with a protective outer container. If the package has been opened outside the sterile field, the product must be considered non-sterile.

## PREPARATION OF ACCESS PORT OR LAP-BAND® SYSTEM

Depending on the area or location of the repair, it is possible that the Access Port and/or the LAP-BAND® System will need to be flushed with sterile saline to remove trapped air prior to placement. Refer to the LAP-BAND® System Surgical Procedures for additional information if necessary.

## RETURNED GOODS POLICY

Authorization must be received from your distributor prior to return of the merchandise. Merchandise returned must have all the manufacturer's seals intact to be eligible for credit or replacement. Products returned may be subject to restocking charges.

## DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Apollo Endosurgery, Inc. product(s) described in this publication. To the fullest extent permitted by applicable law, Apollo Endosurgery, Inc. disclaims all liability for any indirect, special, incidental, or consequential damages, regardless of whether such liability is based on contract, tort, negligence, strict liability, products liability or otherwise. The sole and entire maximum liability of Apollo Endosurgery, Inc., for any reason, and buyer's sole and exclusive remedy for any

cause whatsoever, shall be limited to the amount paid by the customer for the particular items purchased. No person has the authority to bind Apollo Endosurgery, Inc. to any representation or warranty except as specifically set forth herein. Descriptions or specifications in Apollo Endosurgery, Inc. printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties or recommendations for use of the product in specific circumstances. Apollo Endosurgery, Inc. expressly disclaims any and all liability, including all liability for any direct, indirect, special, incidental, or consequential damages, resulting from reuse of the product.

#### **AUTHORIZED TRAINING PROGRAM AND PRODUCT INFORMATION**

LAP-BAND® System placement is an advanced laparoscopic procedure. Surgeons planning LAP-BAND® System placement must participate in a LAP-BAND® System training program authorized by Apollo Endosurgery or an authorized Apollo Endosurgery distributor. This required training program is specific to the Apollo Endosurgery LAP-BAND® System and does not qualify for use with other gastric bands.

**CAUTION:** Federal law (USA) restricts this device to sale by or on the order of a physician.

The LAP-BAND® Adjustable Gastric Banding System contains no latex or natural rubber materials.

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All rights reserved.
worldwide trademarks or registered trademarks of
© 2014 Apollo Endosurgery, Inc., Austin, TX